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Product description:

Manufacturer's information applies to all of Kemmer Nail supplied rotary instruments as carbide / diamond and ceramic cutters used in the nail and foot care and are reusable.

Limited number of reprocessing:

The products described are reusable, but they should not be put into use in wear or damage. Also a frequent reprocessing has absolutely no influence on the service life of the instruments. It is recommended that a check of the instruments before cleaning / sterilization for damage. Under no circumstances instruments should be used if the following damage are recognized:

- tools or instruments with outbreaks at the cut edges
- worn diamond grinder (bare spots)
- instruments with a radial runout mistake

Basically, all instruments must be immediately professionally cleaned / sterilized after use on patients before they can be reused. The storage should be followed in a clean location or in a locked instrument case.

Manual cleaning ultrasonic bath:

Before use of the ultrasonic bath the instruments rinse of dirt under running water and rinse with a suitable cleaning brush. Place instruments then into a suitable container (strainer) and place in the filled with a detergent and disinfectant ultrasonic unit. The cleaning and disinfection must be carried out according to the manufacturer of the cleaning agents and disinfectants. After the reaction time all the instruments are immediately using a full deionized water (DI water) to rinse completely and thoroughly. For drying compressed air is recommended.

After drying, a visual inspection with appropriate magnification should be. Should still be on the instruments residual contamination further mechanical cleaning is required until no more dirt can be seen.

Sterilization:

The described instruments can be done with steam sterilization fractionated vacuum process at 134 °C in a device according to EN285 / EN 13060. To avoid staining and corrosion, the steam must be free of harmful ingredients. The recommended limits for particle contents in feed water and steam condensate are defined by EN 285. The manufacturer's instructions must be observed. Storage of sterilized instruments after steam sterilization should be appropriate under individual packs of dust, moisture and recontamination.

Basic information:

Kemmer Nail guaranteed in accordance with the above procedure described and professionally following manufacturer's instructions of disinfectant and appropriate equipment manufacturers that no damage or shorten the life occurs in the cleaning / sterilizing instruments. All our rotary instruments such as carbide / diamond and ceramic cutters are resistant against all located in the market and approved cleaning agents and disinfectants for such applications as described above.